ARTWORK PENCIL

Millnet Ref: 00580115-001-000
 Job Title: Sinclair Ellansé Face S 2x1ml IFU - Leaflet
 Version: 1 (Packaging Artwork Proof)
 Date: 03-Feb-21

Brand: Sinclair
 Range: Ellansé
 Product: Face
 Variant: S
 Market: Brazil;Korea;Saudi Arabia;Spain

Job Type: IFU - Leaflet
 Pack Size: 2x1ml
 Man. Site: Aqis Medical
 Cutter Ref: AQ-LFT-519.75x210mm

New Item Code:
 Prev Item Code: TR1700004
 Part Number: 1601.168.02
 Prev Part Number: 1601.168.01


Barcode Type:
 Magnification:
 No./NE:

Colours:
 Black

Miller Graphics UK
 2nd Floor, Minster House, 23 Flemingate,
 Beverley, East Yorkshire, HU17 0NP, England
 Tel: +44 (0) 1262 603495
 www.millergraphics.com

UNDER NO CIRCUMSTANCES SHOULD THIS ARTWORK BE ALTERED WITHOUT PRIOR PERMISSION FROM MILLER GRAPHICS

[EN]

DESCRIPTION
 ELLANSÉ is a sterile, latex-free, non-pyrogenic, bioresorbable, non-permanent implant, whose principle component is synthetic Poly-ε-Caprolactone (PCL) microspheres suspended in a gel carrier of phosphate buffered saline (pHEur), hyaluron (pHEur) and carboxymethylcellulose (USP). PCL is a well-known totally bioresorbable soft medical polymer. PCL is used in numerous CE-marked and Food and Drug Administration (FDA) approved commercial bioresorbable product applications for several decades world-wide and has demonstrated an excellent safety profile. ELLANSÉ has a particle size range of 25-50 microns and should be injected with a 27 gauge needle.

INTENDED USE/INDICATIONS
 ELLANSÉ (a Poly-ε-caprolactone-based dermal filler) is an injectable implant, indicated for subdermal implantation in the face for the lasting correction of wrinkles and facial aging signs or conditions. ELLANSÉ uses the body's natural response to encapsulate the microspheres and stimulate the formation of connective tissue. Due to the good biocompatibility of the microspheres, the connective tissue growth, which replaces the gel carrier, is normally completed within a few months. Over time, the imperfection or condition is filled with autologous connective tissue. Due to its unique physical properties the bioresorbability of the PCL microspheres is predictable, controlled and tunable.

MODE OF ACTION
 ELLANSÉ is implanted by subdermal injections with a 27G needle. The homogeneous distribution of the microspheres, its predictable, controlled and tunable bioresorption kinetics, and a specifically designed viscosity which is suitable for subdermal injections, gives ELLANSÉ its physical and mechanical properties and prolonged bioresorbability, which make this implant suitable for filling areas of depressed skin.

CONTRAINDICATIONS
 ELLANSÉ is contraindicated for patients with severe allergies manifested by a history of anaphylaxis, or history or presence of multiple severe allergies. ELLANSÉ is not to be used in patients with known hypersensitivity to any of the components. ELLANSÉ should not be used in the case of acute or chronic skin disease (infection or inflammation). ELLANSÉ should not be used in patients with known susceptibility to keloid formation or hypertrophic scarring. ELLANSÉ should not be administered in case of current cortisone treatment, as the growth of connective tissue fibres might be inhibited. ELLANSÉ should not be injected into the periorbital region (eyelids, under-eye dark circles, crow's feet) or glabella region as there is risk of ocular ischaemic events leading to loss of vision. ELLANSÉ should not be used in patients with active sepsis or infection, active (or history of) autoimmune disease.

WARNINGS
 Do not use if package is opened or damaged. Do not resterilize. Do not use if the syringe end cap or syringe plunger is not in place. ELLANSÉ should only be used subdermally. Avoid superficial injections. Do not over-inject the area to be treated. The safety and effectiveness in the lips has not been established. Special care should be taken to avoid injection into the blood vessels. An introduction into the vasculature may occlude the vessels and could cause infarction or embolism leading to ischemia, necrosis or scarring. This has been reported to occur in the lips, nose, nasal region, nasolabial fold, glabella, forehead, temporal area, cheeks or ocular area. Use of ELLANSÉ at specific sites in which an active inflammatory process (skin eruptions such as cysts, pimples, rashes, or hives) or infection is present should be deferred until the inflammatory process has been controlled. Immediately stop the injection and provide urgent medical care if a patient exhibits any of the following symptoms, including: changes in vision, signs of stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur. ELLANSÉ is packaged for single patient use in one treatment session. Re-use can result in increased risk of infection and decreased effectiveness. Concurrent use of this product in conjunction with other aesthetic procedures has not been established. Physicians injecting Ellansé should have a management strategy in place for the immediate management of a patient that develops ocular pain and impending blindness during injection.

POTENTIAL SIDE EFFECTS / ADVERSE EVENTS
 Injection related reactions, including swelling, redness, pain, itching, discoloration or tenderness may occur at the site of the injection. These usually resolve spontaneously within several days after the injection. As with any implant material possible adverse reactions that may occur include, but are not limited to the following: hypersensitivity, allergic reactions, inflammation, infection, fistula formation, haematoma, seroma, extrusion, induration formation, inadequate healing, skin discoloration, oedema, contusions/bruises, ecchymosis, inadequate or excessive augmentation, loss of correction, interference with local circulation if injected into a blood vessel resulting in vessel laceration, occlusion, embolism, infarction, and abscess at implant site which may result in induration and/or scar formation. Nodules or granulomas may form requiring treatment or removal. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischaemia or cerebral hemorrhage, leading to stroke, skin necrosis, damage to underlying facial structures, nausea and vomiting, lack of ocular movement, significant skin changes, ptosis, exotropia and neurological complications. Adverse events, other than mentioned above, could occur as can be with any medical intervention.

All (potential) adverse side effects and incidents should be reported immediately to the Sinclair Pharma Corporate Safety Group: quality@sinclairpharma.com

PRECAUTIONS FOR USE
 The following precautions should be observed: ELLANSÉ should only be used by health care providers with expertise in the correction of volume defects after fully familiarizing themselves with the product and its complete instruction leaflet. The injection site should be cleaned with an antiseptic and free from inflammation or infection. As with all injections, patients treated with anti-coagulants may run the risk of a haematoma or localized bleeding at the injection site. No studies of interactions of ELLANSÉ with drugs or other substances or implants have been made. As with all transcutaneous procedures, ELLANSÉ injection carries a risk of infection. Standard precautions associated with injectable materials should be followed. Safety of ELLANSÉ for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established. As with all injections, patients who are using medications that can prolong bleeding, such as aspirin or warfarin, may, as with any injection, experience increased bruising, bleeding, or haematoma at the injection site. After use, treatment syringes and needles may be potential biohazards. Handle accordingly and dispose of in accordance with accepted medical practice and applicable local, state and federal requirements. The patient should be informed that he or she should minimize exposure of the treated area to excessive sun and UV lamp exposure and extreme cold weather until any initial swelling and redness has resolved. If laser treatment, chemical peeling or any other procedure based on active dermal response is considered after treatment with ELLANSÉ there is a possible risk of eliciting an inflammatory reaction at the implant site. This also applies if ELLANSÉ is administered before the skin has healed completely after such a procedure. ELLANSÉ is supplied in a syringe ready for use. Never mix ELLANSÉ with other dermal filler or neurotoxin products prior to injection of the device. Use of the supplied needle is recommended. The design, diameter and length have been validated for optimum use with the injection. It is also possible to use injection accessories presented in the table on the end of this leaflet. The maximum recommended lengths of needles and cannulas covering a range of gauge sizes are included. The design, diameter and length have been validated for effective use with the injection. The needle and cannula recommendations are based solely on extrusion testing and clinical judgement should be used to determine the gauge and length most suitable for the application area/depth of treatment. Use of injection equipment other than that recommended in the table of alternatives increases the risk of the accessory and/or the Luer lock fitting becoming detached.

DIRECTIONS FOR USE
 This device should only be used by trained practitioners in the field of soft tissue augmentation. Before treatment with ELLANSÉ, the patient should be informed completely of the indications, contra-indications, warnings, precautions for use, possible side effects and mode of administration of ELLANSÉ. A complete medical history is advised to make sure that the treatment is appropriate. Assess the patient's need for pain management.

Pretreatment photographs are recommended. The following items are required for the percutaneous injection procedure:
 The ELLANSÉ syringe(s)
 27 gauge needle with Luer lock fittings

- Prepare patient for percutaneous injection using standard methods. Local or topical anesthesia at the injection site should be used at the discretion of the physician. The treatment injection site should be marked and prepared with a suitable antiseptic.
- Prepare the syringes of ELLANSÉ and the injection needle(s) before the percutaneous injection. A new injection needle(s) shall be used for each syringe.
- The syringe of ELLANSÉ can be twisted onto the Luer lock fitting of the needle. The needle must be tightened securely to the syringe and primed with ELLANSÉ. If excess ELLANSÉ is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until ELLANSÉ extrudes from the end of the needle. If leakage is noted at the Luer fitting, it may be necessary to tighten the needle, or to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the needle.
- Locate the initial site for the implant. Scar tissue and cartilage may be difficult or impossible to treat. Avoid, if possible, passing through these tissue types when advancing the injection needle.
- The amount injected will vary depending on the site, the extent of the restoration or augmentation desired. ELLANSÉ should be injected subdermally.
- Use a 1:1 correction factor. **Do not overcorrect.**
- Insert needle with bevel down at approximately a 30° angle to the skin. Needle should slide in/under the dermis to the point you wish to begin the injection. This should be easily palpable with the nondominant hand.
- If significant resistance is encountered when pushing the plunger, the injection needle may be moved slightly to allow easier placement of the material or it may be necessary to change the injection needle. Needle jams are more likely with use of needles smaller than 27G.
- Advance the needle in/under the dermis to the starting location. Carefully push the plunger of the syringe to start the injection and slowly inject ELLANSÉ material in linear threads while withdrawing the needle. Continue placing additional lines of material until the desired level of correction is achieved. The needle should be moved backwards (tunneling technique) under pressure, filling the created tunnels. The injection pressure is correct if the implant flows slowly, evenly and without excessive pressure.
- The implant material should be completely surrounded by skin or soft tissue without leaving globular deposits. The injected area may be massaged as needed to achieve even distribution of the implant.

HOW SUPPLIED
 ELLANSÉ is a bioresorbable implant in a ready for use sterile syringe application. ELLANSÉ is supplied in a disposable plastic syringe with sterile content with a Luer-lock tip-cap, single packed in a pouch. Patient record labels are included adhered to each pouch.

SHELF LIFE AND STORAGE
 ELLANSÉ must be used prior to the expiration date printed on the package. Store at ambient temperature 15 – 25°C. Do not freeze, and protect from sunlight. Refrigeration is not required. Do not resterilize ELLANSÉ as this may damage or alter the product. Do not use if the package is damaged.

DISPOSAL
 Used and partially used syringes and injection needles could be bio hazardous and should be handled and disposed of in accordance with facility medical practices and local, state or federal regulations.

WARRANTY
 AQTIS Medical BV warrants that reasonable care has been exercised in the design and manufacture of this product. **THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ITS PARTICULAR PURPOSE.** Handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AQTIS Medical BV's control directly affect the product and the results obtained from its use. AQTIS Medical BV's obligation under this warranty is limited to the replacement of this product and AQTIS Medical BV shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly, arising from the use of this product. AQTIS Medical BV neither assumes, nor authorizes any person to assume for AQTIS Medical BV, any other or additional liability or responsibility in connection with this product.

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